Instructions for use
of **Aniball** & **Aniball INCO**
medical device

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CE

EN CONTENT

1. Intended purpose of Aniball and Aniball INCO	3
2. Contraindications of the use of Aniball and Aniball INCO	3–4
3. Safety precautions for the use of Aniball and Aniball INCO	4–6
4. Aniball	6–11
5. Aniball INCO	11–16

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1. INTENDED PURPOSE OF ANIBALL AND ANIBALL INCO

ANIBALL and ANIBALL INCO are medical devices (hereinafter also referred to as „balloon“) for the training of conscious relaxation and activation (contraction) of the pelvic floor muscles of an adult woman. Both balloons will help to discover the pelvic floor and find out how to work with it effectively in the context of the breath.

Aniball is intended for prenatal preparation, prevention of perineal injuries and stress urinary incontinence.

Aniball INCO is intended for the prevention and treatment of stress urinary incontinence.

Please read the following instructions carefully!

2. CONTRAINDICATIONS OF THE USE OF ANIBALL AND ANIBALL INCO

The use of the Aniball/Aniball INCO balloon is limited by the woman's health. Consult a physician, physiotherapist, or midwife to assess the suitability of the exercise and to consider contraindications:

Do not use Aniball/Aniball INCO in the following cases:

- high-risk pregnancy
- any vaginal inflammation (infectious or non-infectious)
- any vaginal bleeding
- presence of condyloma or other infectious pathology in the vagina or vaginal entrance area
- vaginal injury
- any bladder inflammation
- uterine prolapse
- genital herpes
- precancers and tumour pathologies of the cervix, vagina or external genitalia

EN

DE

FR

IT

ES

EN

- during the postnatal puerperium and within the first six weeks following gynaecological interventions and operations (including genital laser treatment)

DE

FR

In the following cases, you need to be extra cautious and consult your doctor:

IT

ES

- presence of varicose veins in the vagina and external genitalia
- low-lying placenta
- reduced sensitivity of the vagina or external genitalia (neurological disease, painkillers administration, alcohol and drug abuse)
- after treated vaginal inflammation (vaginal fragility persists)
- after vaginal surgery (scarring may negatively affect the vaginal elasticity)
- lichen sclerosus of the external genitalia (chronic skin disease)
- cancer in the pelvic area

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user is established.

3. SAFETY PRECAUTIONS FOR THE USE OF ANIBALL AND ANIBALL INCO

- Before the first, and after each subsequent use, wash the silicone balloon thoroughly using oil-free pH-neutral soap and warm water - at least 40°C / 104°F - for 1–2 minutes and rinse under running water. Subsequently, it is advisable to disinfect the balloon using a product suitable for medical silicone. Perform the disinfection according to the manufacturer's instructions. Or boil the balloon in boiling water for at least 5 minutes. Then, allow it to dry on a clean and smooth surface. You may then place the balloon onto the plastic adapter. The balloon fits to the adapter relatively tightly, but it can be disconnected and reconnected at any time.

- After each use, the balloon must be thoroughly washed with pH neutral oil-free soap, rinsed under warm running water, disinfected with a suitable product and stored in a (the original) bag.
- If pauses between exercise sessions exceed 5 days, wash and disinfect or boil the balloon again before the next use.
- If in doubt about choosing a suitable soap and disinfectant, visit our website: **www.aniball.uk / www.aniball.co.nz / www.aniball.com.au**
- Before inserting the balloon into the vagina, it is necessary to check its integrity and correct functionality by inflating it using the pump to 20 pumps. **Perform this check before each use! Never use the product in the case of any damage or signs of material degradation! In case of damage to the balloon, a spare part can be purchased.**
- For comfortable and safe insertion, it is advisable to lubricate the surface of the balloon **using exclusively water-based lubricants suitable for internal use (contact with mucous membranes). Do not use products or condoms containing oils.** Non-water-based lubricants may damage the silicone, causing the balloon to burst. Please be aware that such damage will not be considered a product defect in the case of a warranty claim.
- If you practice any kind of pelvic massage before exercising your pelvic floor muscles with Aniball, use a suitable water-based lubricant instead of massage oil.
- Aniball is recommended for use **from the end of the 36th week** of pregnancy until delivery.
- **Do not try to work through the pain nor intentionally induce pain while exercising. Pain is not a sign of correct exercise. What you should feel is a pleasant pull from the stretching or the muscle activation.**
- **It is not necessary to increase the balloon size at all costs during prenatal preparation with Aniball.** Clinical data shows that for effective prevention of birth injuries, it is sufficient to reach a balloon circumference of **25 cm (approximately 10 pumps)**. **We recommend that you do not go beyond this to avoid injury!**
- **The balloon should never slip out of the vagina too quickly, so hold it by the plastic neck.** A quick slip indicates a loss of control and could result in injury.
- If you feel a sharp or disproportionate pressure and tension in your lower abdomen, take a break or stop exercising.
- During exercise, it is a common mistake **to pull the belly button** into the abdomen or **to pull the gluteus muscles** towards you. This is an inappropriate substitute that prevents the true deep activation of your pelvic muscles.

EN

DE

FR

IT

ES

- Prevent contact between the pump and water when washing the balloon, as water could prevent its proper functioning.
- Always let the balloon dry thoroughly after each wash.
- Store the device separately in the original bag in a dry place and do not expose it to radiant heat or sunlight!
- Prevent the device from making contact with sharp objects.
- Do not use in water.
- Keep the device out of reach of children.
- For safety and hygienic reasons, the device is intended for use by one person only!

If you have any questions, problems or suggestions, please do not hesitate to contact us at: hello@aniball.uk / info@aniball.co.nz / info@aniball.com.au. We are always here for you and ready to help.

4. ANIBALL

Intended purpose


The Aniball medical device is intended for prenatal preparation, prevention of perineal injuries and stress urinary incontinence.

Properties

Aniball is a proven medical device that effectively helps to practice the conscious relaxation of the pelvic floor in the context of breathing without forced pushing, allows you to test the appropriate birth position and pressure, positively stimulates the vaginal muscles and improves perineum elasticity during prenatal preparation. The result is that the expulsion phase of the newborn is facilitated, thus helping to prevent the perineum from tearing or cutting. The use of the device after the sixth week postpartum contributes to the strengthening of the pelvic floor muscles, which leads to a reduction in incontinence in the later period after childbirth.

Aniball consists of:

1. Inflatable balloon made of medical grade silicone
2. Mechanical pump
3. Tube
4. Connecting adapter
5. Two-way valve

I. ANTENATAL PREPARATION WITH ANIBALL

Pre-exercise Information

Aniball is suitable to use from the **end of the 36th week** of pregnancy. Exercise for approximately 15–30 minutes daily. This time and frequency of exercise are only indicative, each woman requires individual time for proper relaxation. At the beginning of the exercise, **take the position that suits you best**; for example, lying down on your back with slightly raised knees, kneeling, lying on your side, standing or crouching. Childbirth in the position in which your body feels most comfortable is usually the safest.

Inserting the Balloon

Settle into your chosen position and give yourself a few moments to become calm and relaxed. During exercise, remind yourself to remain relaxed and avoid muscle tension or muscle spasms in your neck, arms, hands or buttocks. The ability to stay calm and relaxed (not only during the following exercise but also during labour) will have a significant positive impact on your birth experience.

The silicone balloon can be inflated for 1–2 pumps to give it some resistance. It is also possible to insert the balloon completely deflated. To insert comfortably and safely, it is recommended to moisten the surface of the balloon using exclusively a water-based lubricant intended for internal use (see **Safety Precautions for Use**). Now take hold of the balloon behind the plastic neck and gently push the larger rounded part into the vagina.


EN

DE

FR


IT

ES

Once inserted, if the balloon does not cause any pain or unpleasant feeling, you can begin to inflate it. A pump with a valve at the end of the tubing is provided for inflation. During inflation, the volume of the balloon increases and the pressure inside the vagina increases with it. This pressure is caused by the inflation of the balloon against the vaginal wall, stretching the surrounding muscles. You will experience very similar pressure during your childbirth.

How much inflation is needed is individual to your body. The first time you use Aniball, it should be to a **maximum of five pumps**, including the pumps you performed before insertion. Subsequently, as you become comfortable with Aniball, you can increase the number of pumps, but **always with great care and always stop if you feel pain. Whenever the balloon seems too big or unpleasant, you can use the valve to release air and thus reduce its volume.**

Turn the small wheel lightly using two fingers. Stop turning as soon as you encounter resistance. Do not tighten it. The valve is extremely sensitive and works as a seal.

Exercise Instructions

Remain in the position that you selected for the exercise.

1. Conscious Control Over the Movement of the Balloon

You will learn to control the balloon's movements inside your vagina. This requires some conscious effort and focus on your body and it takes practice to feel confident. The balloon

can, depending on the activity of your pelvic floor muscles, move within the vagina, either out or inwards. If the pelvic floor muscles are actively withdrawn, the balloon sits deeper, but if you manage to relax these muscles consciously, the balloon will have the opportunity to move closer to the vaginal entrance. Training the controlled movement out of the vagina will help to prepare for your upcoming delivery.

You can now try these conscious movements. While exhaling, try to squeeze the balloon using your pelvic floor muscles and feel its movement (now it should be submerging deeper into the vagina). Then inhale and relax your muscles (feel the balloon moving out of the vagina). The important thing is to gain control over these movements and to realize that only pelvic muscles allow the balloon (as well as the baby) to move down the birth canal outwards. The aim of this exercise is not to expel the balloon fully out of the vagina yet. Repeat the entire exercise of pulling the balloon in and out ten times.

2. Expulsion of the Balloon

Relax again and gradually move on to the expulsion of the inserted balloon.

Use what you have learned in the previous exercise. While inhaling, relax your pelvic floor muscles and feel the balloon coming out of the body **without any forced pushing**. **The aim of the exercise is to expel the balloon using only your breathing (inhaling) and conscious relaxation of the pelvic muscles.**

When exhaling, simply try to keep the balloon in the same position as upon the inhaling phase, and do not draw it deeper into the vagina. With another inhale, release your pelvic floor and move the balloon one step lower. You are now consciously expelling the balloon until the perineum slightly stretches out and the balloon comes out of the vagina. **The balloon should not slip out of the vagina too quickly; therefore, hold it by the plastic neck** (just as you do during insertion). If the balloon slips out too fast, it indicates a loss of control and could lead to injury. Even your baby will be held during delivery, so feel free to hold Aniball.

The pressure, which you consciously create throughout the expulsion, **should only come from the pelvic floor area**. Tensioning of the neck, shoulders and hands, as well as contracting the abdomen and contracting the buttocks, is not desirable and indicates an incorrect function of the pelvic floor during the expulsion. If your body automatically does this, don't worry, it's very good that you've started practising with Aniball to gently prepare for birth.

EN

3. Expansion of the Balloon

DE

You have started exercising after the end of the 36th week of pregnancy, and you are approaching the date of your delivery with each passing day. It is therefore advisable to gradually get used to expelling a larger circumference of the balloon.

FR

However, be careful! It is not beneficial to increase the size of the balloon if it causes pain or significant discomfort. Clinical data shows that for effective prevention of birth injuries, it is sufficient to reach a balloon circumference **of 25 cm (approximately 10 pumps).** **We recommend not to go beyond this to avoid injuries!!** If you are worried about the balloon being too large or you do not feel good about it and you are comfortable with the current smaller size, it is perfectly all right, and there is no need to force yourself into increasing balloon size. The most important thing is **conscious control over the movements of the balloon inside the vagina and getting used to the birth canal for the child's future passage.** The same applies when the expulsion of a larger balloon causes great pain or even the possibility of injury to surrounding tissue. In this case, there is also no reason to force yourself into increasing balloon size.

IT

ES

Follow the rule to always exercise with such a size that is tolerable and relatively pleasant both for your body and your mind. Never take on more than you can handle, do not expose your body to pain or injury. **Do not compare yourself with other mums** at the same stage of pregnancy. Each of us is unique and it's important to listen to, and respect your body's limitations. Remember that the hormones that are released within your body during child-birth will have a huge impact on your body's responses - so not being able to practice with a larger balloon size during pregnancy does not mean you will struggle to deliver your baby. You may use the enclosed tape measure to keep track of the size of the balloon. But once again, **numbers are not your goal! The goal is to be a healthy, calm and ready mother who can control the pelvic floor to encourage an easier passage for your baby with minimal injury to your own body.**

II. USE ANIBALL TO EXERCISE AFTER THE SIXTH SUNDAY AFTER BIRTH TO STRENGTHEN THE PELVIC FLOOR TO PREVENT INCONTINENCE

Every woman should know that after a vaginal birth and a caesarean section, it is necessary to continue working with the pelvic floor. Not only childbirth, but the whole pregnancy can

cause a disruption of the integrity of the pelvic floor, which leads to serious problems with its function.

Consistent use of the balloon after the sixth week after pregnancy positively contributes to restoring the proper function of the pelvic floor muscles and preventing incontinence. The classic **Aniball** or better, its smaller variant **Aniball INCO**, can be used for this purpose. **Exercising with Aniball after the sixth week postpartum is the same as exercising with Aniball INCO to prevent incontinence. Instructions for use and exercise procedure after the sixth week can be found below.**

5. ANIBALL INCO

Intended purpose


Aniball INCO is intended for the prevention and treatment of stress urinary incontinence.

Properties

Aniball INCO (hereinafter also referred to as “balloon”) is used to train the proper function (activation/relaxation) of the pelvic floor muscles for the prevention and treatment of stress incontinence. The balloon inserted into the vagina is easily inflated with air to the desired comfortable size. Thanks to its variable shape and size, it adapts optimally to the individual anatomy of each woman. The balloon inserted in this way must not cause pain or other discomfort. In case the balloon pressure is perceived unpleasantly or painfully, it is possible to reduce it at any time by releasing the air valve on the pump body. However, light balloon pressure is desirable to feel its position and resistance during exercise. By activating and, conversely, releasing the pelvic floor muscles, there is a slight movement of the balloon inside the vagina (inwards of the vagina and outwards), which provides very important feedback on the quality of work of the pelvic floor muscles.

Aniball INCO consists of:

1. Inflatable balloon
made of medical grade silicone
2. Mechanical pump
3. Tube
4. Connecting adapter
5. Two-way valve

EN

Safety information

DE

FR

IT

ES

- Women who have given birth or are suffering from incontinence should consult their gynaecologist, urologist or physiotherapist regarding the individual suitability of the exercise.
- If you have decided to use **Aniball** or **Aniball INCO** after giving birth, keep in mind that you can **exercise at the earliest after the end of the sixth week and any birth injuries have healed.**
- **During exercise, never hold your breath! (Except before the intentional coughing during the pressure resistance training).**
- Sometimes, you may feel fatigue in your pelvic muscles. That's perfectly all right; it indicates you are exercising correctly. If the fatigue is so great that the muscles cannot be activated, stop the exercise.
- **Never allow an inserted and simultaneously inflated balloon to be expelled directly out of the vagina during exercise. We always recommend that you release it to the minimum size using the valve before removing it.** Along with the breath, relax the pelvic floor muscles, grasp the balloon by the plastic neck and slowly pull it out of the vagina, or slide it outwards so that it slides out spontaneously.

Recommendations for Using Aniball INCO

- In order to use this device effectively and reap the rewards of your efforts, you need to exercise regularly and with persistence. Regular exercise (every day for 5–15 minutes) alleviates the difficulties within approximately three months.
- Learn to control your pelvic floor muscles during the day, even when not exercising with the balloon. Before coughing or sneezing, imagine you are squeezing the balloon and preventing it from sliding out of the vagina. This activation will help you to regulate urine leakage in a strenuous situation. **Never “strengthen” your pelvic floor by interrupting urination! This incorrect procedure increases the risk of urinary tract infection, and may lead to disruption of the micturition reflex which contributes to incontinence!**

Exercise Instructions

You should exercise at least once a day for 5–15 minutes. All exercises should be performed calmly, slowly and as closely to these instructions as possible. Exercise your pelvic floor muscles with an empty bladder.


1. Body Position During Exercise

Start with a position lying down on your back, bend your knees and anchor your feet on the ground. Place your feet at hipwidth distance apart. If you feel your head is leaning back, place a pillow or a folded towel under it, so the cervical spine is straight, without a backward bend. Your back should rest comfortably on the mat without a bend or arch in your lower back. Leave your hands loosely alongside your body or place them comfortably on your lower abdomen.

2. Inserting the Balloon

Before the balloon is inserted, relax and do not rush. The vagina may be dry and sensitive; therefore, we recommend **applying appropriate lubricant** (see **Safety Precautions for Use**). You may leave the balloon deflated or slightly inflate it for 1–2 pumps. Then, grab the balloon behind the plastic neck and insert it into the vagina. If you have inserted the balloon inflated and now feel uncomfortable pressure inside your vagina, you can reduce its size using the valve at the end of the pump to achieve a comfortable volume. If you have inserted the balloon but do not feel its slight resistance, we recommend inflating it to a size which is comfortable for you but also makes you feel slight resistance of the balloon against the vaginal walls. Once you reach optimal feeling without uncomfortable pressure or pain, but at the same time a slight resistance of the balloon, you can start with the exercise. Perform the exercises below with respect to your own body.

If you are tired or have a headache, take a break at any time during the exercise or stop with the exercise for that day. You must not feel pain during any of these exercises.

EN

3. Breathing Exercise

DE

In this exercise, you will learn to activate and release the pelvic floor muscles in coordination with your breathing. First, slow down your breath and be aware of individual inhales and exhales.

FR

IT

ES

While inhaling, imagine your breath flowing through your chest down to your abdomen and lower abdomen. Visualize the inhaled air reaching the pelvic floor muscles, all the way to the place of the inserted balloon. If your palms are placed on your abdomen or groin, you can feel the growing inner volume in this area which aims to gently push out the balloon.

While **inhaling**, teach your pelvic floor muscles to consciously relax in order to feel the balloon moving under the gentle pressure in the outward direction. However, you should not loosen your muscles too much to avoid expelling the balloon completely out of the vagina (the goal is a conscious and controlled release).

Upon each **exhale**, learn to firmly grip the balloon and move it upward deeper inside the vagina (a conscious activation of the pelvic muscles). At the same time, the internal volume of the abdomen created by the previous inhale may shrink again. It is very useful to visualize the balloon squeezed by the vaginal walls, not by the external urethral or rectal sphincters. Try not to squeeze your buttock muscles when trying to hug the balloon, nor to draw your belly button into your belly. **While exhaling, leave the abdominal area relaxed without clenching it** (learn to use only your pelvic muscles).

Inhale again with conscious relaxing of muscles and release of the balloon, similarly to the previous exercise. According to your own natural breathing cycle, alternate between inhale (releasing the grip and gently pushing the balloon downward) and exhale (active grip of the balloon and slight upward movement). The duration of this exercise depends on your concentration and the fatigue rate of your muscles. Initially, try to repeat each inhale-exhale cycle at least five times; gradually, you may train yourself up to ten repetitions.

4. Prolonging the Grip Time

During this exercise, you will increase the strength and endurance of your pelvic muscles.

Perform the exercise in the same way as the previous one (see **Breathing Exercise**) with one difference in the activation stage (grip of the balloon), where you consciously prolong the grip of the balloon to hold it for at least ten seconds (you may count to ten). While holding the balloon, try either to do long and continuous exhales, or to inhale again and

exhale (but, if possible, try it without losing the grip of the balloon). Repeat the entire cycle of continuous activation 5–10 times.

5. Resistance Against Increased Pressure

The aim of this exercise is to manage situations of excessive abdominal pressure, which can lead to unwanted release of the pelvic floor or urine leakage.

Take the position of lying on your back as in previous exercises. While inhaling, lead your breath to the lower abdomen and groin to increase the volume of this area. Now try squeezing the balloon upon inhaling. Exhale and try to not lose the obtained volume and consciously keep the balloon squeezed using your pelvic muscles. Repeat 5–10 times.

Continue with exhaling while trying to cough at the same time, without letting the balloon slide down or all the way out of the vagina. Repeat 3–5 times.

6. Relaxation

*Get ready for the most popular and calming part of your exercise dedicated to relaxing and loosening your pelvic muscles after their previous strengthening. You can deflate the balloon first and remove it carefully from the vagina (see **Removing the Balloon**). Alternatively, you can leave the balloon inside the vagina; this is not wrong. With the balloon still inserted, you will feel the relaxing of your pelvic muscles more intensely; however, do not let the inserted and inflated balloon slip all the way out of the vagina during the exercise.*

You should remove the balloon, or let it slip out freely, always while deflated!

Stay on your back, pull your knees up to your chest with your hands, so your feet no longer touch the floor, and assume the position of a “curled ball”. Let your head rest on a thin pillow or folded towel. Now focus on the inhaling phase and imagine the entire area of your pelvic floor expanding and extending in all directions, like a blooming flower. While exhaling, stop squeezing the balloon, and simply relax. Repeat the inhale and let your pelvic floor once again open up like a flower with no signs of a grip or forced squeezing. Dedicate a lot of time for the relaxing exercise in order to give your muscles enough time to relax. Therefore, repeat the entire inhale-exhale cycle 15–30 times.

EN

7. Removing the Balloon

DE

FR

IT

ES

Before removing the balloon from the vagina, we recommend deflating it to a minimum size using the valve. Inhale and relax your pelvic floor muscles, then take the balloon by the plastic neck and slowly pull it out of the vagina. Alternatively, slide it outward, so it slips out on its own. Continue with washing the balloon, see **Safety Precautions for Use**.

You can find more information at: www.aniball.uk
www.aniball.co.nz / www.aniball.com.au

Shelf life: indicated on the medical device label and packaging.

Prior to disposal, clean and disinfect in accordance with the instructions for use!

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